

JET

JET

uprighter JET,



5- 10-

| | |
|---|---|
| 855-512 - Uprighter Bayonet - Left 10 | |
| 855-511 - Uprighter Bayonet - Right 10 | |
| Button Attachment 852-170 qty 10 | |
| 855-606 240 gm NITI Springe qty 10 | 855-606 180 gm NITI Springe qty 10 |
| 855-608 Activation Lock - qty 10 | 855-510 Distal Stop - qty 10 |
| .050" Hax Wrench - qty 1 855-515 | |
| 855-602 - Bayonet Director - qty 10 | |




Dr. Aldo n and Sig. Mauro Testa

the uprighter **JET**

AMERICAN ORTHODONTICS

UPRIGHTER JET

855-513



①

②

③

④

⑤

⑥

(.1).

() (.2).

(.).



.040" .045" (1.15)




, 240 ,


180 .

| | 180 | .012 | .055 | | | | |
|-----|-----|------|------|----|-----|-----|-----|
| () | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
| Min | 20 | 40 | 62 | 76 | 105 | 126 | 180 |
| | 17 | 33 | 52 | 69 | 94 | 115 | |


| | 240 | .014 | .055 | | | | |
|-----|-----|------|------|-----|-----|-----|-----|
| () | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
| Min | 32 | 62 | 101 | 140 | 172 | 205 | 240 |
| | 21 | 55 | 90 | 125 | 158 | 195 | |



.1



.2



.3



①

②

③

④

4

2-

⑤